

Millipore®

Filtration, Separation
& Preparation

Samplicity® G2 Filtration System for HPLC Samples

The better way to use Millex® filters



The life science business of Merck KGaA,
Darmstadt, Germany operates as
MilliporeSigma in the U.S. and Canada.


**MILLIPORE
SIGMA**



HPLC sample filtration is here to stay.

According to a 2015 survey by LC/GC magazine, sample filtration was performed by a larger percentage of separation scientists than extraction, centrifugation, precipitation, or concentration*.

With downstream detection methods becoming more and more sensitive to particulate contamination, filtration is as crucial as ever.



*Based on survey of 413 respondents. Each respondent was asked to choose multiple methods that they were currently using.

Love Millex® filters, but have tired thumbs?

For chromatographers who filter dozens of samples a day through Millex® syringe filters, especially for hard-to-filter, viscous or particulate-laden samples, traditional filtration can lead to manual fatigue, and even repetitive stress injury. Nevertheless, they stick with Millex® filters, which are specified in numerous analytical methods and standard operating procedures.

Introducing the Samplicity® G2 Filtration System for Millex® Filters.

By adding a simple adapter funnel to standard Millex® 33 mm filters, we've enabled easy, vacuum-driven filtration of up to 8 samples at once, directly into HPLC vials, in seconds.

Look up the Millex® 33 mm filters you use for manual syringe filtration in the table below to find the corresponding Samplicity® G2 filter catalog numbers.



Millex® Filter Description	Catalog No.* (Manual)	Catalog No. (Samplicity® G2 system)
Millex® Filters for the Samplicity® G2 System, 0.22 µm Durapore® PVDF	SLGV033N-	SAMP2GVNB
Millex® Filters for the Samplicity® G2 System, 0.45 µm Durapore® PVDF	SLHV033N-	SAMP2HNVB
Millex® Filters for the Samplicity® G2 System, 0.22 µm Millipore Express® PLUS PES	SLGP033N-	SAMP2GPNB
Millex® Filters for the Samplicity® G2 System, 0.45 µm Millipore Express® PLUS PES	SLHP033N-	SAMP2HPNB
Millex® Filters for the Samplicity® G2 System, 0.20 µm Nylon	SLGN033N-	SAMP2GNNB
Millex® Filters for the Samplicity® G2 System, 0.45 µm Nylon	SLHN033N-	SAMP2HNNB

*The last character of the Millex® part number represented by the "-" may be an S, K, or B.

80%

of chromatographers report using MilliporeSigma filters for sample filtration

Enjoy the Quality of Millex® Filters.

You trust Millex® filters for HPLC sample prep because of our highest quality design, manufacturing, materials and QC standards.

Millex® filters are known for their:

- Minimal hold-up volume
- Filter and seal integrity
- Low extractables, low binding


Millex® filters are available in a variety of membranes and polypropylene housings to ensure chemical compatibility with a range of samples and solvents.

The adapter funnels provided for the Smplicity® G2 system are made of the same polypropylene as our Millex® filter housings. Together, the filter plus adapter are HPLC-certified to be low in extractable impurities. That means you'll continue to enjoy clear chromatography results, without the pain of syringe filtration.

Experience the efficient way to use Millex® filters: vacuum filtration.

Any process that uses Millex® 33 mm filters can now be performed using the vacuum-driven, time-saving, Smplicity® G2 filtration system – even for small-volume (as little as 300 µL) or hard-to-filter samples.

If manual filtration is used, hard-to-filter samples can require excessive manual force, greater time per sample, and can yield lower recovery – all of which may reduce laboratory efficiency. In contrast, the Smplicity® G2 filtration system filters even highly viscous samples in seconds, with minimal manual effort.



Hard to filter? No More! Efficient filtration of viscous samples using the Smplicity® filtration system with 0.45 µm PVDF Millex Smplicity® Filters. Filtration times for 1 mL of indicated samples are shown. Results represent the average of four filtrations, and error bars represent standard deviation.



Save your thumbs.

According to a Baseline Risk Identification of Ergonomic Factors (BRIEF™) survey, the Smplicity® G2 system reduced the ergonomic impact of filtration compared to manual syringe filtration, as shown in the chart below.

Smplicity® G2 System has lower ergonomic impact than manual syringe filtration

Sample Type	Job Hazard Score		
	Low (0-9)	Medium (10-29)	High (30-49)
Syringe Filtration of Non-viscous Sample (water)		12.8	
Syringe Filtration of Viscous Sample (16% ketchup)		16.0	
Smplicity® G2 System Filtration of Viscous Sample (16% ketchup)	8.0		

*For reference, a task involving manipulating a heavy (13 kg) tool was found to have a medium/high job hazard score of 27. Note that the job hazard score for Smplicity® G2 remains at 8 regardless of the viscosity of the sample.



Ultra-convenience

In addition to our Millex® 33 mm syringe filters, Millex Smplicity® filters for the Smplicity® G2 system take convenience in sample filtration one step further. These all-in-one filters feature:

- Integrated funnels for easiest sample loading
- Filters in perforated strips of 4 for fastest setup
- Quality membrane options, including PVDF (0.45 µm), hydrophilic PTFE (0.45 µm) and hydrophilic PTFE (0.2 µm)



See it in action: Smplicity® G2 Videos

Watch a quick video how-to, plus a side-by-side comparison with syringe filtration.

Just visit:
www.emdmillipore.com/smplicity

Ordering information for the Smplicity® G2 System

Millex® Filter Description	Quantity	Catalog No. (Smplicity® G2 system)
Millex® Filters for the Smplicity® G2 System, 0.22 µm Durapore® PVDF	250	SAMP2GVNB
Millex® Filters for the Smplicity® G2 System, 0.45 µm Durapore® PVDF	250	SAMP2HVNB
Millex® Filters for the Smplicity® G2 System, 0.22 µm Millipore Express® PLUS PES	250	SAMP2GPNB
Millex® Filters for the Smplicity® G2 System, 0.45 µm Millipore Express® PLUS PES	250	SAMP2HPNB
Millex® Filters for the Smplicity® G2 System, 0.20 µm Nylon	250	SAMP2GNNB
Millex® Filters for the Smplicity® G2 System, 0.45 µm Nylon	250	SAMP2HNNB
Smplicity® G2 Filtration System, Bold Blue	1	SAMP2SYSB
Smplicity® G2 Filtration System, Glossy Green	1	SAMP2SYSG
Millex Smplicity® Filters, 0.45 µm hydrophilic PTFE	96	SAMPLCR01
Millex Smplicity® Filters, 0.45 µm hydrophilic PTFE	384	SAMPLCR04
Millex Smplicity® Filters, 0.20 µm hydrophilic PTFE	96	SAMPLG001
Millex Smplicity® Filters, 0.20 µm hydrophilic PTFE	384	SAMPLG004
Millex Smplicity® Filters, 0.45 µm hydrophilic PVDF	96	SAMPHV001
Millex Smplicity® Filters, 0.45 µm hydrophilic PVDF	384	SAMPHV004
Smplicity® G2 Lid	1 pk	SAMP2LID
Chemical Duty Pump, 115V/60Hz	1	WP6111560
Chemical Duty Pump, 220V/50Hz	1	WP6122050
Chemical Duty Pump, 100V/50-60Hz	1	WP6110060

*Note: Millex® filters are available in packs of 1000. Please visit: www.emdmillipore.com

Upgrade your current Smplicity® filtration system to the Smplicity® G2 model!

Get all the support you need for the HPLC workflow.

- Sample preparation
- Columns
- Mobile phases
- Technical information

Visit: www.emdmillipore.com/hplc-workflow



To place an order or receive technical assistance

In the U.S. and Canada, call toll-free 1-800-645-5476

For other countries across Europe and the world,
please visit: www.emdmillipore.com/offices

For Technical Service, please visit:
www.emdmillipore.com/techservice

MilliporeSigma.com

