




Mycotoxins

- Mycotoxins are ISO 17034 accredited
- Mix includes 5 of the most commonly tested mycotoxins
- Custom standards are available



🌐 EmeraldScientific.com
✉ Sales@EmeraldScientific.com
☎ (877) 567-3598 x1

Spex CertiPrep is an
Antylia Scientific company.
Find out more at antylia.com.



Mycotoxins

Mold and fungi are all over the world and are of major concern in agricultural products. Fungi and their spores can spread throughout crops and food stores tainting crops. Fungi can proliferate into colonies with the right environmental conditions which destroy crops by consuming them or poison them with toxic secondary metabolites called mycotoxins.

Spex CertiPrep offers an analytical standard mix of five of the most commonly tested mycotoxins (four aflatoxins and ochratoxin A) for your mycotoxin testing needs. Our mycotoxins are ISO 17034 accredited. Custom standards are also available. Contact the Emerald Scientific team at +1.7877.567.3598 x1 or via email at Sales@EmeraldScientific.com.

Mycotoxins Standard Mix					
Compound	CAS #	Concentration	Matrix	Volume	Part #
Aflatoxin B1	1162-65-8	10 µg/mL of each compound	Acetonitrile	1 mL	MYCO-1
Aflatoxin B2	7220-81-7				
Aflatoxin G1	1165-39-5				
Aflatoxin G2	7241-98-7				
Ochratoxin A	303-47-9				



Organic Certified
Reference Materials



Supplied with a
Certificate of Analysis



ISO Accredited
Standards



🌐 EmeraldScientific.com
✉ Sales@EmeraldScientific.com
☎ (877) 567-3598 x1

Spex CertiPrep is an
Antylia Scientific company.
Find out more at antylia.com.

