

CleanPrep™ Dead Air Box

110V



GENERAL SPECIFICATIONS		
Construction Material	Work Surface	Polypropylene
	Shell	Polypropylene
	Front Sash	Polycarbonate
	Pipette Shelf	Polycarbonate
Controls		ON/OFF Switch
Lighting		Fluorescent
UV Light		254-nm Shortwave UV Light
UV Light Timer		0-15 Minute UV Light Timer

DIMENSIONS

Product #	External Dimensions			Internal Dimensions			Sash Opening Dimensions		
	Width	Depth	Height	Width	Depth	Height	Width	Opening 1	Opening 2
MY-DB24	24.0"	24.0"	29.5"	23.5"	20.0"	29.0"	22.0"	10.5"	23"
MY-DB36	36.0"	24.0"	29.5"	35.5"	20.0"	29.0"	34.0"	10.5"	23"
MY-DB48	48.0"	24.0"	29.5"	47.5"	20.0"	29.0"	46.0"	10.5"	23"

TECHNICAL SPECIFICATIONS

Product #	Electrical			Unit Weight*
	Voltage	Amps	Power	
MY-DB24	110V	0.5 A	50 W	60 lbs
MY-DB36	110V	0.5 A	50 W	75 lbs
MY-DB48	110V	0.7 A	80 W	90 lbs

*Unit weight is an approximation



MYSTAIRE®

Science.
Safety.
Solutions.

www.mystaire.com • contact@mystaire.com

Phone: (919) 229-8511 • Toll Free: (877) 328-3912 • Fax: (919) 528-1744