

MY-PCR Prep Station

110V



GENERAL SPECIFICATIONS		
Construction Material	Work Surface	Polypropylene
	Shell	Polycarbonate
	Front Sash	Polycarbonate
	Filtration Housing	ABS
	Pipette Shelf	Polycarbonate
Controls		ON/OFF Switch
Blowers		Brushless-Sparkless
Air Quality		Class 100 / ISO 5
Lighting		Fluorescent
UV Light		254-nm Shortwave UV Light
UV Light Timer		15-Minute UV Light Timer
HEPA Filter Alarm		Audible / Visual

DIMENSIONS

Product #	External Dimensions			Internal Dimensions			Sash Opening Dimensions		
	Width	Depth	Height	Width	Depth	Height	Width	Opening 1	Opening 2
MY-PCR24	24"	24"	30"	23"	23"	18"	21"	10"	20"
MY-PCR32	32"	24"	30"	31"	23"	18"	29"	10"	20"
MY-PCR48	48"	24"	30"	47"	23"	18"	45"	10"	20"

TECHNICAL SPECIFICATIONS

Product #	Electrical			Pre-filters Required	Main Filters Required	Unit Weight*
	Voltage	Amps	Power			
MY-PCR24	110V	1.2 A	130 W	1	1	55 lbs
MY-PCR32	110V	1.2 A	130 W	1	1	70 lbs
MY-PCR48	110V	2.1 A	230 W	2	2	95 lbs

\*Unit weight is an approximation.



MYSTAIRE®

Science.  
Safety.  
Solutions.

[www.mystaire.com](http://www.mystaire.com) • [contact@mystaire.com](mailto:contact@mystaire.com)

Phone: (919) 229-8511 • Toll Free: (877) 328-3912 • Fax: (919) 528-1744