

MY-ISO5 Vertical Laminar Flow Cabinet

110V



GENERAL SPECIFICATIONS		
Construction Material	Work Surface	Polypropylene
	Shell	Polycarbonate
	Front Sash	Polycarbonate
	Filtration Housing	ABS
Controls		ON/OFF Switch
Blowers		Brushless-Sparkless
Air Quality		Class 100 / ISO 5
Lighting		Fluorescent
IV Bar		Optional, Stainless Steel
HEPA Filter Alarm		Audible / Visual

DIMENSIONS

Product #	External Dimensions			Internal Dimensions			Sash Opening Dimensions		
	Width	Depth	Height	Width	Depth	Height	Width	Opening 1	Opening 2
MY-VLF24	24"	24"	30"	23"	23"	18"	21"	10"	20"
MY-VLF32	32"	24"	30"	31"	23"	18"	29"	10"	20"
MY-VLF32T	32"	24"	39"	31"	23"	27"	29"	12"	
MY-VLF48	48"	24"	30"	47"	23"	18"	45"	10"	20"
MY-VLF48T	48"	24"	42"	47"	23"	29"	45"	10"	

TECHNICAL SPECIFICATIONS

Product #	Electrical			Pre-filters Required	Main Filters Required	Unit Weight*
	Voltage	Amps	Power			
MY-VLF24	110V	1.0 A	110 W	1	1	55 lbs
MY-VLF32	110V	1.0 A	110 W	1	1	70 lbs
MY-VLF32T	110V	1.0 A	110 W	1	1	80 lbs
MY-VLF48	110V	1.9 A	200 W	2	2	95 lbs
MY-VLF48T	110V	1.9 A	200 W	2	2	105 lbs

*Unit weight is an approximation.



MYSTAIRE®

Science.
Safety.
Solutions.

www.mystaire.com • contact@mystaire.com

Phone: (919) 229-8511 • Toll Free: (877) 328-3912 • Fax: (919) 528-1744