

Simplify your centrifuge selection process with the OHAUS guide to centrifuge and rotor bundles! Featuring combinations designed around specific Research or Clinical applications, OHAUS has the right centrifuge and rotor packages to help you **Do More** in your laboratory.

Research Applications

High-Speed Micro Centrifuge Bundles



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
1	High-Speed Micro, Non-Refrigerated, Standard Rotor	5515 Centrifuge (30130867)	1	15,200	21,953	24
		30130870 (24 × 1.5/2.0ml Rotor)	1			



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
2	High-Speed Micro, Refrigerated, Standard Rotor	5515R Centrifuge (30130869)	1	15,200	21,953	24
		30130870 (24 × 1.5/2.0ml Rotor)	1			



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
3	High-Speed Micro, Non-Refrigerated, Biocontainment Rotor	5515 Centrifuge (30130867)	1	13,500	19,150	30
		30130872 (30 × 1.5/2.0ml Rotor)	1			



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
4	High-Speed Micro, Refrigerated, Biocontainment Rotor	5515R Centrifuge (30130869)	1	13,500	19,150	30
		30130872 (30 × 1.5/2.0ml Rotor)	1			

High-Speed Micro Centrifuge Bundles



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
5	High-Speed Micro, Non-Refrigerated, High Volume	5515 Centrifuge (30130867)	1	15,000	21,379	44
		30372717 (44 × 1.5/2.0ml Rotor)	1			



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
6	High-Speed Micro, Refrigerated, High Volume	5515R Centrifuge (30130869)	1	15,000	21,379	44
		30372717 (44 × 1.5/2.0ml Rotor)	1			

Microplate Bundle



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
7	Microplate Bundle	FC5714 Centrifuge (30314811)	1	4,500	2,716	6 plates
		30314824 (2 × 3 Microplate Rotor)	1			

Low Volume—Cell Culture Bundles



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
8	Low Volume Cell Culture Bundle, Non-Refrigerated	FC5706 Centrifuge (30332132)	1	6,000	4,427	6
	Fixed Angle	30130878 6 × 50mL Rotor (included with above item#)	0			
		30130893 (15mL Adapter), Pk/2	3			



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor	
9	Low Volume Cell Culture Bundle, Non-Refrigerated	FC5718 Centrifuge (30314813)	1	5,000	3,885		
	Swing-Out	30314823 (4 × 100ml Swing-Out Rotor)	1				
		30314882 (3 × 15mL bucket)	1				6
		30314881 (2 × 50mL bucket)	1				4



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor	
10	Low Volume Cell Culture Bundle, Refrigerated	FC5718R Centrifuge (30314815)	1	5,000	3,885		
	Swing-Out	30314823 (4 × 100ml Swing-Out Rotor)	1				
		30314882 (3 × 15mL bucket)	1				6
		30314881 (2 × 50mL bucket)	1				4

Mid-Volume—Cell Culture Bundles



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
11	Mid-Volume Cell Culture Bundle, Non-Refrigerated	FC5816 Centrifuge (30314817)	1	4,500	3,780	
	Swing-Out	30314820 (4 × 250mL Swing-Out Rotor)	1			
		30314913 (7 × 15mL bucket)	1			14
		30314906 (3 × 50mL bucket)	1			6



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
12	Mid-Volume Cell Culture Bundle, Refrigerated	FC5816R Centrifuge (30314819)	1	4,500	3,780	
	Swing-Out	30314820 (4 × 250mL Swing-Out Rotor)	1			
		30314913 (7 × 15mL bucket)	1			14
		30314906 (3 × 50mL bucket)	1			6

High Volume—Cell Culture Bundles



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
13	High Volume Cell Culture Bundle, Non-Refrigerated	FC5916 Centrifuge (30553037)	1	4,000	3,434	
	Swing-Out	30553085 (4 × 750mL Swing-Out Rotor)	1			
		30553121 (Bucket 1 × 750ml D99mm FB)	2			
		30553130 (Rack 7 × 50ml D29mm FA)	1			
		30553134 (Rack 17 × 15ml D17mm FA)	1			34



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
14	High Volume Cell Culture Bundle, Refrigerated	FC5916R Centrifuge (30553103)	1	4,500	4,346	
	Swing-Out	30553085 (4 × 750mL Swing-Out Rotor)	1			
		30553121 (Bucket 1 × 750ml D99mm FB)	2			
		30553130 (Rack 7 × 50ml D29mm FA)	1			
		30553134 (Rack 17 × 15ml D17mm FA)	1			34

Clinical Applications

Low Volume—Clinical Bundles



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
1	Low Volume 5/7mL Clinical Bundle, Non-Refrigerated	FC5707 Centrifuge (30393189)	1	6,800	4,445	8
	Fixed Angle	30472305 8 × 15mL fixed angle rotor (included with above item#)	0			
		30130889 Adapter, D13.5mm, 7ml Tubes × 2	3			
		or				
		30130890 Adapter, D13.5mm, 5ml Tubes × 2	3			



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
2	Low Volume 5/7mL Clinical Bundle, Non-Refrigerated	FC5706 Centrifuge (30130876)	1	6,000	4,427	12
	Fixed Angle	30130877 (12 × 15mL fixed angle rotor)	1			
		30130889 Adapter, D13.5mm, 7ml Tubes × 2	3			
		or				
		30130890 Adapter, D13.5mm, 5ml Tubes × 2	3			

Low Volume—Clinical Bundles



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
3	Low Volume 5/7mL Clinical Bundle, Non-Refrigerated	FC5714 Centrifuge (30314811)	1	4,500	2,830	30
	Fixed Angle	30314832 (Rotor, Angle, 30 × 15ml RB/FA, ID)	1			
		303130889 (Adapter, D13.5mm, 7ml Tubes × 2)	15			
		or				
		303130890 (Adapter, D13.5mm, 5ml Tubes × 2)	15			



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
4	Low Volume 5/7mL Clinical Bundle, Non-Refrigerated	FC5714 Centrifuge (30314811)	1	4,500	3,350	28
	Swing-Out	30314822 (4 × 100mL swing out rotor)	1			
		30314857 (Rack, 7 × 5/7mL Vacutainer), 2/pk	2			



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
5	Low Volume 10mL Clinical Bundle, Non-Refrigerated	FC5714 Centrifuge (30314811)	1	4,500	3,350	20
	Swing-Out	30314822 (4 × 100mL swing out rotor)	1			
		30314858 (Rack, 5 × 10mL Vacutainer), 2/pk	2			

Mid Volume—Clinical Bundles



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
6	Mid-Volume 5/7mL Clinical Bundle, Refrigerated	FC5718R Centrifuge (30314815)	1	4,500	3,350	28
	Swing-Out	30314822 (4 × 100mL swing out rotor)	1			
		30314857 (Rack, 7 × 5/7mL Vacutainer), 2/pk	2			



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
7	Mid Volume 10mL Clinical Bundle, Refrigerated	FC5718R Centrifuge (30314815)	1	4,500	3,350	20
	Swing-Out	30314822 (4 × 100mL swing out rotor)	1			
		30314858 (Rack, 5 × 10mL Vacutainer), 2/pk	2			



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
8	Mid Volume 3/5/7mL Clinical Bundle, Non-Refrigerated	FC5816 Centrifuge (30314813)	1	4,500	3,780	40
	Swing-Out	30314820 (4 × 250mL swing out rotor)	1			
		30314917 (Rack, 10 × 1.6-7mL), 2/pk	2			

Mid Volume—Clinical Bundles



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
9	Mid Volume 3/5/7mL Clinical Bundle, Refrigerated	FC5816R Centrifuge (30314819)	1	4,500	3,780	40
	Swing-Out	30314820 (4 × 250mL swing out rotor)	1			
		30314917 (Rack, 10 × 1.6-7mL), 2/pk	2			



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
10	Mid Volume 10mL Clinical Bundle, Non-Refrigerated	FC5816 Centrifuge (30314817)	1	4,500	3,780	32
	Swing-Out	30314820 (4 × 250mL swing out rotor)	1			
		30304367 (Rack, 8 × 10mL), 2/pk	2			



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
11	Mid Volume 10mL Clinical Bundle, Refrigerated	FC5816R Centrifuge (30314819)	1	4,500	3,780	32
	Swing-Out	30314820 (4 × 250mL swing out rotor)	1			
		30304367 (Rack, 8 × 10mL), 2/pk	2			

High Volume—Clinical Bundles



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
12	High Volume 3/5/7mL Clinical Bundle, Non-Refrigerated	FC5916 Centrifuge (30553037)	1	4,000	3,434	108
	Swing-Out	30553085 (4 × 750mL swing out rotor)	1			
		30553121 (Bucket 1 × 750ml D99mm FB)	2			
		30553139 (Rack 27 × 1.8/7ml D13mm Vac/Sar × 2)	2			



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
13	High Volume 3/5/7mL Clinical Bundle, Refrigerated	FC5916R Centrifuge (30553103)	1	4,500	4,346	108
	Swing-Out	30553085 (4 × 750mL swing out rotor)	1			
		30553121 (Bucket 1 × 750ml D99mm FB)	2			
		30553139 (Rack 27 × 1.8/7ml D13mm Vac/Sar × 2)	2			

High Volume—Clinical Bundles



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
14	High Volume 10mL Clinical Bundle, Non-Refrigerated	FC5916 Centrifuge (30553037)	1	4,000	3,434	72
	Swing-Out	30553085 (4 × 750mL swing out rotor)	1			
		30553121 (Bucket 1 × 750ml D99mm FB)	2			
		30553140 (Rack 18 × 4/10ml D16mm Vac/Sar)	2			



Bundle #	Description	Contents	QTY	RPM	RCF	Tubes per Rotor
15	High Volume 10mL Clinical Bundle, Refrigerated	FC5916R Centrifuge (30553103)	1	4,500	4,346	72
	Swing-Out	30553085 (4 × 750mL swing out rotor)	1			
		30553121 (Bucket 1 × 750ml D99mm FB)	2			
		30553140 (Rack 18 × 4/10ml D16mm Vac/Sar)	2			